



Superior Class B & S Autoclaves

THE PROMEDCO RANGE OF AUTOCLAVES ARE OF THE HIGHEST STANDARD, BACKED BY A GENEROUS TWO YEAR WARRANTY FOR TROUBLE-FREE OPERATION.

PROMEDCO IS A WHOLLY OWNED AUSTRALIAN COMPANY SERVICING ALL CAPITAL CITIES AND REGIONAL AREAS.

THE PROMEDCO SUPERIOR RANGE OF AUTOCLAVES ARE DESIGNED FOR TROUBLE-FREE OPERATION. COMPETITIVELY PRICED WITH GUARANTEED RELIABILITY, THE PROMEDCO RANGE ARE FAST BECOMING THE PREFERRED CHOICE OF DENTISTS, DOCTORS, PODIATRISTS, BEAUTICIANS, TATTOOISTS AND VETS.

SIMPLE & EASY TO USE WITH ALL THE FUNCTIONALITY YOU NEED

Our autoclaves are user friendly and meet the most stringent requirements for hollow and wrapped loads. You can choose between a quick start and a pre-set, custom-defined option. With an inbuilt printer and USB

WHY WE'RE DIFFERENT

Your new autoclave will come fully validated at no extra charge. In addition we provide you with a starter kit of steri pouches, class 5 & 6 indicators and a sterilisation log book. After your purchase we can provide all your consumables at extremely competitive prices.

FEATURES AT A GLANCE



Fully Compliant

ProMedCo guarantees all our autoclaves are compliant.
Your safety is our priority.



100% Satisfaction

We are so confident in our products, we proudly offer a 30 day money back guarantee. If you are not 100% happy with the unit, we will organise for it to be collected and a full refund to be provided.



Training & Support

We offer in-house or skype training sessions on your new autoclave, along with full telephone support.



Service

We service & validate most major brands of autoclaves. Annual validation is a legal requirement for autoclaves used in skin penetration industries. If your Autoclave requires service and can't be repaired onsite we will provide you with a loan unit to allow your business to continue uninterrupted.



BEST IN CLASS AUTOCLAVES



This sleek design comes with a door that covers access to the printer roll and water outlets.



One piece solid stainless steel chamber with three trays and easy access to water tank for filling and cleaning purposes.



THOMAS PUMPS: GERMAN BUILT AND DESIGNED TO LAST



The ULKA Water Pump

Our autoclaves contain the industry gold standard of Thomas pumps. Thomas pumps are synonymous with high quality, durability and performance. Thomas pumps feature numerous filters, high quality valves and ULKA water pumps, ensuring that your autoclave won't miss a beat.



	8L ProMedCo Superior Autoclave	12L ProMedCo Superior Autoclave	18L ProMedCo Superior Autoclave	22L ProMedCo Superior Autoclave
Class	IIb	IIb	IIb	IIb
Certification	CE 0197	CE 0197	CE 1097	CE 0197
Power Requirement	10A, 240V AC	10A, 240V AC	10A, 240V AC	10A, 240V AC
Wattage	2.4 KW	2.4 KW	2.4 KW	2.4 KW
Sterilization temperature	121 - 134 °C	121 - 134 °C	121 - 134 °C	121 - 134 °C
Max working pressure	0.21-0.23MPA	0.21-0.23MPA	0.21-0.23MPA	0.21-0.23MPA
Diameter of chamber	17CM	20CM	25CM	25CM
Depth of chamber	32CM	36CM	35CM	45CM
Volume of chamber	8L	12L	18L	23L
Material of chamber	Stainless steel 304	Stainless steel 304	Stainless steel 304	Stainless steel 304
Colours available	White	White	White	White
Screen type	LCD	LCD	LCD	LCD
Screen size	6 x 4.25cm	6 x 4.25cm	6 x 4.25cm	6 x 4.25cm
Language available	English, Italian, German, Chinese	English, Italian, German, Chinese	English, Italian, German, Chinese	English, Italian, German, Chinese
Lock system of door	Mechanical & electronic lock	Mechanical & electronic lock	Mechanical & electronic lock	Mechanical & electronic lock
User password	Yes	Yes	Yes	Yes
Cycle Recording	Thermo printer & USB	Thermo printer & USB	Thermo printer & USB	Thermo printer & USB
Net weight	35kg	51kg	53kg	55kg
Shipping weight /carton	42kg	59kg	61kg	63kg
Product size	56 x 44 x 36cm	61 x 50 x 43cm	61 x 50 x 43cm	65 x 50 x 43cm
Shipping crate size	65 x 45 x 42cm	73 x 57 x 51cm	73 x 57 x 51cm	77 x 57 x 51cm



CALL US FOR A FREE CONSULTATION

02 9734 0673

AUTOCLAVE | VALIDATION | DELIVERY | CONSUMABLES